# miniland

## thermotalk plus

Español	3	Português	35	Deutsch	64	Polski	94
English	19	Français	49	Italiano	79	Русский	109

## ÍNDICE

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## INDEX

1. Introduction
2. Contents
3. Safety instructions
4. Quick start guide
5. Product features
6. Instructions for use
7. Advanced settings
8. Product care
9. Troubleshooting
10. Specifications & symbols

## 1. INTRODUCTION

Congratulations on your purchase of this clinical IR thermometer with pulse reading function.

Before using this device, please read all instructions in this manual carefully. For the safe and correct use of this product, please read the safety precautions explained in this user manual.

Please store for future reference as it contains important information.

## 2. CONTENTS

- 1 clinical IR thermometer
- 1 stand
- 1 storage and travel bag
- 2 x AAA batteries
- 1 instruction manual and guarantee

If any of the above is missing, please contact your retailer.

## 3. SAFETY INSTRUCTIONS

### 3.1. WARNING

Warning marks and symbols are included for safety and a correct use of this product and prevention of any injury.

See below for the description of warning marks and symbols:

- **WARNING:** Means the possibility of personal injury in case of improper use.
- **NOTICE:** Means the possibility of injuries or damage\* in the case of improper use.

[\*This covers any damage to your house, personal belongings and pets]

- **PROHIBITION:** Means that the detailed items expressed in words or figures within or beside the symbol are totally forbidden. The symbol on the left means forbidden in general.

- **MUST BE OBSERVED:** This means obligatory with the detailed items expressed in words or figures within or beside the mark. The one on the left means that it is compulsory in general.

- **THE MEANING OF THE SYMBOL:** Means Type-BF device.

• The patient may be an operator and can change the batteries under normal circumstance and maintain the device and its accessories according to the user manual.

• It is dangerous for patients to perform a self-evaluation and self-treatment based on the results of the reading. Be sure to follow your doctor's instructions.

*\*A self-evaluation may cause deterioration in the state of the patient's health.*

• **Don't touch or blow on the infra-red sensor.**

*\*A dirty infra-red sensor may cause an inaccurate reading.*

• **Clean the infra-red sensor with a soft wet cloth.**

*\* Do not clean with toilet paper or a tissue as this may scratch the infrared sensor and cause an inaccurate reading.*

• **Clean the infra-red sensor every 4 uses (following the instructions given on the manual, section 8.1 CARE AND MAINTENANCE).**

• **Keep the machine out of children's reach**

*\* If a child takes his own reading, it may cause damage to his ear. Please consult the doctor at once if a battery is accidentally swallowed.*

• **When the device is used to measure the temperature of a child, it has to be operated by an adult. Adults can measure their own temperature.**

• **In the case of having performed any outdoor activities, no reading should be taken until the patient's body temperature is the same as the room temperature.**

*\*If the patient has used a water-filled pillow or hot water bottle, this may cause an inaccurate reading.*

• **If there is a difference in temperature between the storage area and the place where the reading is to take place, leave the thermometer at room temperature for approximately 30 minutes.**

*\*A failure to do so may cause an inaccurate reading.*

• **If the reading exceeds 38°C and / or the patient feels uncomfortable, the user should get in contact with a doctor.**

• **Stop using the product in the event of any pain.**

*\*It is likely to cause damage to the inside of the patient's ear.*

• **It is not recommended to use the product by users suffering ear disease including otitis and tympanitis.**

*\*The affected part is likely to get worse.*

• **Don't use the product in a wet ear.**

*\*It may cause damage to the ear.*

• **Do not dispose of the batteries by throwing them into a fire.**

*\* They may explode.*

**FEVER!** If the reading exceeds 38°C (100.4°F) and / or the patient feels uncomfortable, the user should get in contact with a doctor.

• It is dangerous for patients to perform a self-evaluation and self-treatment based on the results of the reading. Be sure to follow your doctor's instructions.

• Keep the thermometer out of children's reach.

- The device is intended for occasional reading and monitoring of body temperature in a domestic environment. Do not use the product for other purposes.

### 3.2. RECOMMENDATIONS

- The product should not be exposed to any chemical solvent, direct sunshine or high temperatures.
- Do not drop, step on or subject the product to any vibration or impact.
- Do not use a mobile phone near the thermometer when it is being used

### IMPORTANT

- Changes and modifications carried out by anyone other than the MINILAND technical service will not be covered by the product guarantee.
- The warranty excludes misuse:
  - Unauthorized disassembly and modification of the product.
  - If the product is dropped during application or transportation.
  - An incorrect use of the product.

## 4. QUICK START GUIDE

### 4.1. TAKING A FOREHEAD TEMPERATURE

• **Ensure the forehead probe is fitted to the unit**, then locate the probe gently in the middle of the forehead. Press the forehead button [] once. The device will beep twice to indicate that the reading is complete (1-2 seconds). The LCD will display the temperature. You will also hear the result via the speaker if the TALK function is enabled.

• **FEVER!** If the reading exceeds 38°C (100.4°F) and / or the patient feels uncomfortable, the user should get in contact with a doctor.

• You can press the heart rate button [] at any time and this will return you to the time mode.

### 4.2. TAKING IN AN INNER EAR TEMPERATURE

• **Remove the forehead probe cover from the unit**, and gently insert the temperature probe as deep as possible into the ear. Press the ear button [] once. The device will beep twice to indicate that the reading is complete. The LCD will display 'ooo' and the temperature. You will also hear the result via the speaker if the TALK function is enabled.

• **FEVER!** If the reading exceeds 38°C (100.4°F) and / or the patient feels uncomfortable, the user should get in contact with a doctor.

• When the "ooo" disappears from the screen, the next reading can be taken.


• A warning sound will be heard if you press any button before all the "ooo" disappear.

## 5. PRODUCT FEATURES

### Intended for use:

The device is intended for the occasional reading and monitoring of human body temperature in a domestic environment.


### 5.1 PARTS

### IMPORTANT

1. Ear button
2. Probe
3. Stand
4. LCD display
5. Heart rate button
6. Talk button
7. Forehead probe cover
8. Forehead button
9. Battery compartment
10. Infrared sensor


## 5.2. LCD DISPLAY

1. AM/PM
2. Time, temperature or pulse display
3. Memory number
4. Memory number or year display
5. Month display
6. Date display
7. Heart rate reading symbol
8. Ear temperature reading symbol
9. Speaker
10. Forehead temperature reading symbol
11. Set display
12. Body temperature reading symbol
13. Low battery symbol
14. Pulse
15. Cleaning indicator

## 5.3. FEATURE DESCRIPTION

### 5.3.1. ROOM TEMPERATURE, DATE & TIME DISPLAY

The LCD display will alternate every 5 seconds between showing the time and then the room temperature.

The date range is 2001-2099 and the time is displayed in a 24 hour clock format.

### 5.3.2. READING AND DISPLAY OF ROOM TEMPERATURE

Reading ranges from 5.0°C to 59.9°C (41.0°F to 139.8°F).

### 5.3.3. READING AND DISPLAY OF BODY TEMPERATURE

Readings range from 32.0°C to 42.9°C (89.6°F to 109.3°F).

### 5.3.4. MEMORY

Maximum of 15 sets of temperature records can be stored in the devices memory. The temperature is stored together with date, memory number and mode (Ear or Forehead).

### 5.3.5. AUDIO

Temperature, time and records can be heard using the TALK function.

### 5.3.6. HEART RATE MONITOR

Range: 40~140bpm.

## 6. INSTRUCTIONS FOR USE

### 6.1. LANGUAGE SETUP

To change between languages, press and hold the talk button [TALK] for 3 seconds and release. Repeat to scroll through Spanish, English, German, French, Italian, Russian or mute function. When the TALK function is enabled (i.e. you have selected a language), the [TALK] icon is displayed.

### 6.2. PREPARATIONS

1. For a correct reading, the thermometer must be in the same environment 30 minutes before the reading takes place.

Replace any low-voltage battery to ensure full power supply.

2. When using for the first time, we recommend taking the reading 3 times in the same ear. The maximum reading should be recorded in the case of different results as there may be a slight difference of between (+/- 0.3°C, +/- 0.5°F) if the thermometer is used continuously.

3. Any sudden changes in the temperature of the environment are likely to affect an accurate reading.
4. Three readings can be taken consecutively. Wait only 30 seconds before taking another reading on the same patient to avoid excessive cooling of the skin. However, after the 3rd reading, wait 10 minutes before taking another reading.
5. For users in sound health, both ears are almost the same in terms of body temperature. You should select the ear that shows a constantly higher temperature (feels warmer).
6. The reading must be taken in a clean ear to assure maximum accuracy. Please clean the ear before taking a reading.
7. The reading must be taken when the patient is in a state of rest for maximum accuracy.
8. If the patient has used a water-filled pillow or hot water bottle, this may cause an inaccurate reading.
9. The temperature is likely to rise immediately after exercise or bathing. In the case of having performed any outdoor activities, no reading should be taken until the patient's body temperature is the same as the room temperature. Allow 30 minutes for the body to cool down before taking a temperature.
10. If there is a difference in temperature between the storage area and the place where the reading is to take place, leave the thermometer at room temperature for approximately 30 minutes, failure to do so may cause an inaccurate reading.

**NOTE:** The unit may not function correctly in an environment with electrostatic discharge so it may be necessary to reset the unit.

## 6.3. TAKING A TEMPERATURE

### 6.3.1. TAKING A FOREHEAD TEMPERATURE

1. **Ensure the forehead probe is fitted to the unit,** then locate the probe gently in the middle of the forehead. Press the forehead button [  ] once. The device will beep twice to indicate that the reading is complete (1-2 seconds). The LCD will display the temperature. You will also hear the result via the speaker if the TALK function is enabled.

2. **FEVER!** If the reading exceeds 38°C (100.4°F) and / or the patient feels uncomfortable, the user should get in contact with a doctor.

- The result will automatically be recorded in the memory. Once 15 recordings have been stored, on taking the 16th recording, the first set of data will be removed and the new recording will take the top position.
- Three readings can be taken consecutively. However, you should wait for at least 10 minutes before taking a new reading.
- Display of the results:
  - An accurate value display for result ranges between 32.0°C and 42.9°C (89.6°F and 109.3°F)
  - "HI" is displayed for result above 42.9°C~109.3°F
  - "LO" is displayed for result below 32.0°C~89.6°F
- Once the reading has finished, the unit will automatically return to time mode after 30 seconds.
- The display will automatically turn off after 3 minutes after the last key has been pressed (if sleep mode is enabled).
- You can press the heart rate button [♥] at any time and this will return you to the time mode.

### 6.3.2. TAKING IN AN INNER EAR TEMPERATURE

1. **Remove the forehead probe cover from the unit**, and gently insert the temperature probe as deep as possible into the ear. Press the ear button [👂] once. The device will beep twice to indicate that the reading is complete. The LCD will display '000' and the temperature. You will also hear the result via the speaker if the TALK function is enabled.

2. **FEVER!** If the reading exceeds 38°C (100.4°F) and / or the patient feels uncomfortable, the user should get in contact with a doctor.


- The result will automatically be recorded in the memory. Once 15 recordings have been stored, on taking the 16th recording, the first set of data will be removed and the new recording will take the top position.
- Three readings can be taken consecutively. However, you should wait for at least 10 minutes before taking a new reading.
- Display of the results:
  - An accurate value display for result ranges between 32.0°C and 42.9°C (89.6°F and 109.3°F)
  - "HI" is displayed for result above 42.9°C~109.3°F
  - "LO" is displayed for result below 32.0°C~89.6°F
- Once the reading has finished, the unit will automatically return to time mode after 30 seconds.
- The display will automatically turn off after 3 minutes after the last key has been pressed (if sleep mode is enabled).

- When the “ooo” disappears from the screen, the next reading can be taken.
- A warning sound will be heard if you press any button before all the “ooo” disappear.

#### 6.4. HEART RATE MONITOR

- To monitor your heart rate, press the heart rate button [♥]. The device will beep once and the LCD will display the heart symbols. Lightly rest your thumb or index finger on the heart rate button (do not depress the button). When the device detects the heart rate, the symbol [♥] will beat to the beat of your pulse. After approx 15 seconds you will hear a beep to indicate that the reading is complete. Your pulse rate will be displayed on the LCD.
- If the heart symbol [♥] does not flash, please make sure your finger or thumb is placed correctly on the heart button and that you are pressing down hard enough. **NOTE:** Excessive pressure will obstruct the local blood circulation and the device will not be able to detect the heart rate.

- Press the forehead button [⌘] to exit heart rate mode at any time.

#### Caution:


- Do not keep the [♥] pressed down when the device is taking a reading of your heart rate.
- If the symbol [♥] does not flash, please make sure your finger is placed correctly and you are pressing down hard enough. Excessive pressure will obstruct the local blood circulation and the device will not be able to detect the heart rate.
- The maximum reading time is 15 seconds.

- Avoid swimming or intense physical exercise at least half an hour before measuring heart rate.
- Do not measure heart rate if your hand is shaking.
- You must put your finger on the heart rate button [♥] without moving it.
- Whilst measuring, do not talk or move.
- If the measurement is incorrect, you may take it again. However, you must wait 5 to 15 minutes before taking the next measurement.
- Make sure you use the same finger if you decide to repeat the heart rate measurement so that you can compare the results.
- The results obtained after measuring heart rate must only be regarded as a reference. They cannot be used as the basis for carrying out clinical diagnoses or evaluations. It is dangerous to draw conclusions or start treatments based on these results. Please consult your doctor and follow his indications if you don't feel well.

## 6.5. MEMORY MODE

Each temperature result will automatically be stored in the memory. Once 15 temperatures have been stored, on taking the 16th temperature, the first result will be removed and the new recording will take the top position.

- Press the forehead [⌘] and talk button [🗨️] together once to enter memory mode. "--" will be displayed if no data is stored.
- To scroll between stored data, press the forehead button [⌘] whilst in memory mode. Press and hold the forehead button [⌘] to scroll quickly between stored data.
- Press the talk button [🗨️] to hear the currently selected stored temperature data. Press any button to stop the sound during the announcement.
- Maximum storage capacity is 15 sets of data with date and memory number shown.
- The unit automatically returns to clock mode 5 seconds after the last button has been pressed.
- You can press the heart rate button [♥] at any time before the 30 seconds is up and this will return you to the time mode.

## 7. ADVANCED SETTINGS

### 7.1. FAHRENHEIT/CELSIUS SETUP

To change between Fahrenheit (°F) and Celsius (°C) press the ear button [👂] and talk button [🗨️] together.

CELSIUS

FAHRENHEIT

### 7.2. TIME & DATE SETUP

1. Press and hold the forehead [👁️] and talk [🗨️] buttons together for approximately 3 seconds to enter the setup mode. The hour digits will flash indicating you are in setup mode.

2. The device will automatically exit setup mode 30 seconds after the last button has been pressed.

3. Press the heart rate button [❤️] at any time to exit setup.

4. Press the talk button [🗨️] to select between setting the HOUR, MINUTE, YEAR, MONTH, DAY and SLEEP functions.

HOUR

MINUTE

YEAR

MONTH

DAY

5. Press forehead button [👁️] to modify the selected function. Pressing and holding the forehead button [👁️] will scroll through the digits more quickly.

6. Press the talk button [🗨️] to confirm the settings and go to the next function.

### 7.3. SLEEP MODE SETUP

It is recommended that you use sleep mode as this will prolong the battery life. In sleep mode the display will automatically turn off after 3 minutes to save battery power if sleep mode is set to "ON".

1. Press and hold the forehead [  ] and talk [  ] buttons together for approximately 3 seconds to enter the setup mode. The hour digits will flash indicating you are in setup mode.

2. The device will automatically exit setup mode 30 seconds after the last button has been pressed.

3. Press the heart rate button [  ] at any time to exit setup.

4. Press the talk button [  ] to cycle through the HOUR, MINUTE, YEAR, MONTH and DAY functions till you reach the SLEEP function.

SLEEP ON

SLEEP OFF

5. Press forehead button [  ] to modify the selected function.

6. Press the talk button [  ] to confirm the settings and exit setup.

## 8. PRODUCT CARE

### 8.1. CARE AND MAINTENANCE


1. Clean the thermometer and replace the cover to protect the lens after each use.
2. Remove any stains on the casing with a soft dry cloth.
3. The device is not waterproof - do not submerge in water. It should only be cleaned with a damp cloth.
4. Do not use detergents containing abrasive elements or benzene.
5. Clean the thermometer and replace the cover to protect the lens after each use.
6. Do not expose the thermometer to direct sunshine, high temperatures, damp conditions, fire, flame, vibration or impact.
7. Remove the battery for long time of non use.
8. Keep the probe surface clean, otherwise the reading will not be accurate.
9. Do not touch or blow on the infrared sensor as a dirty infrared sensor may cause an inaccurate reading.
10. Do not clean the infrared sensor with toilet paper or a tissue as this may scratch it and cause an inaccurate reading.
11. Clean the probe with clean, soft, paper moistened in water or alcohol, and take the reading after the water on the surface of the probe all evaporates.

## 8.2. WARNING OF DIRT ON THE SENSOR

In order to remind the user to clean the infrared sensor, this thermometer has a special alarm. If the thermometer hasn't been used for a long time or if several measures have been taken, the LCD will display the cleaning symbol to remind the user to clean the infrared sensor. The thermometer will return back to standby mode after pressing any button.

## 8.3. CALIBRATION

- This thermometer is calibrated at the time of manufacture. If the thermometer is used according to instructions, periodic recalibration is not required. If at any time you question the accuracy of the measurement, please contact the retailer immediately.
- The unit may not function correctly in an environment with electrostatic discharge so it may be necessary to reset the unit.
- Don't attempt to modify or reassemble the thermometer.

## 8.4. BATTERIES

### 8.4.1. REPLACE THE BATTERIES

- Replace batteries if low battery symbol is displayed [  ].
- Remove the battery compartment flap.
- Use only new batteries (2 x AAA battery).
- Replace the battery compartment firmly.
- Remove the batteries if storing the product for long periods of non use.

### 8.4.2. BATTERY & PRODUCT DISPOSAL INFORMATION

- Remember to dispose of the product responsibly.
- Do not mix products bearing the crossed-out bin symbol with your general household waste.
- For the correct collection and treatment of these products take them to the collection points designated by your local authority. Alternatively contact the retailer who sold you the product.
- Please consult the doctor at once if a battery is accidentally swallowed.
- Do not dispose of the batteries by throwing them into a fire as they may explode.
- For the disposal of batteries and waste, please act according to local law.

## 9. TROUBLESHOOTING

PROBLEMS	CHECK LIST	SOLUTIONS
No response/ Automatic reset when the insulator is taken out.	Battery used up.	Insert a new battery.
	Battery in wrong polarity.	Remove battery, turn it over and re-insert.
	Poor battery contact.	Remove battery, clean contacts and re-insert.
Battery symbol  on LCD.	Low battery life.	Insert a new battery.
“Lo” on screen (result below 32.0°C or 89.6°F).	Ensure thermometer is placed correctly in the eardrum / forehead.	Follow user’s manual to point at the eardrum / forehead.
“Hi” on screen (result below 42.9°C or 109.3°F).	Please check the operation method.	Follow user’s manual for correct reading.
“Errp” on screen.	Hardware problem.	Consult your dealer
“ErrE” on screen.	Environment temperature beyond the scope of : 10°C~40°C (50.0°F~104.0°F) for ear. 15°C~40°C (59.0°F~104.0°F) for forehead.	Condition the thermometer until it reaches environment temperature: 10°C~40°C (50.0°F~104.0°F) for ear. 15°C~40°C (59.0°F~104.0°F) for forehead.
“Err” on screen.	Operation error, more than 15 seconds without test result in heart rate mode.	Please repeat the operation, following the instructions in user manual.

## 10. SPECIFICATIONS & SYMBOLS

### 10.1 TECHNICAL SPECIFICATION

- Product Name: Clinical IR thermometer
- Model: 89068
- Power supply: DC3V (2 x AAA battery)
- Mode of operation: Adjusted mode
- Reference body site: Core
- Temperature sensing Part: IR sensor
- Range: 32.0°C ~ 42.9°C (89.6°F ~ 109.3°F)
- Room Temperature: 5.0°C ~ 59.9°C (41.0°F ~ 139.8°F)

#### Laboratory accuracy:

- +/- 0.2°C (or +/- 0.4°F)
- 35.5°C ~ 42.0°C (or 95.9°F ~107.6°F)
- +/- 0.3°C (+/- 0.5°F)
- 35.0°C ~ 35.4°C (or 89.6°F ~ 95.8°F)
- 42.1°C ~ 42.9°C (or 107.7°F ~ 109.3°F)

- Room Temperature: +/- 2°C (4°F)
- Size: 127(L) x 39(W) x 63(D) mm
- Weight : 67g (without batteries)

### Working environment:

10°C~40°C (50.0°F~104.0°F) for Ear  
 15°C~40°C (59.0°F ~104.0°F) for Forehead

Relative humidity: 15%RH~85%RH

• Transportation and storage environment: -20°C~55°C (-4°F~131°F)

Relative humidity: 15%RH~90%RH

Atmospheric pressure: 700hPa-1060hPa

• Precision/repeatability: +/- 0.3°C (+/- 0.5°F)

## 10.2. NORMALIZED SYMBOLS

The following is a list of symbols that can be found on the product, packaging and instructions.

- Attention: Read the operating Instructions for use!

- Applied part of type BF

- Disposal in accordance with Directive 2012/19/UE (WEEE)

- Complies with the European Medical Device Directive (93/42/EEC) and amended by directive 2007/47/EC requirements. Notified body TÜV Rheinland (CE0197)

- Serial number

- The name and the address of the manufacturer.

- Authorized Representative in the European Community

### IMPORTANT INFORMATION REGARDING ELECTRO MAGNETIC COMPATIBILITY (EMC)

With the increased number of electronic devices such as PC's and mobile (cellular) telephones, medical devices in use may be susceptible to electromagnetic interference from other devices. Electromagnetic interference may result in incorrect operation of the medical device and create a potentially unsafe situation.

Medical devices should also not interfere with other devices.

In order to regulate the requirements for EMC (Electro Magnetic Compatibility) with the aim to prevent unsafe product situations, the EN60601-1-2 standard has been implemented. This standard defines the levels of immunity to electromagnetic interferences as well as maximum levels of electromagnetic emissions for medical devices.

This medical device conforms to this EN60601-1-2 standard for both immunity and emissions. Nevertheless, special precautions need to be observed:

- The device needs to be installed and put into service in accordance with the information provided in the user manual.
- The use of accessories and cables other than those specified by Manufacturer, with the exception of cables sold by Manufacturer as replacement parts for internal components, may result in increased emission or decreased immunity of the device.
- The medical devices should not be used adjacent to or stacked with other equipment. In case adjacent or stacked use is necessary, the medical device should be observed to verify normal operation in the configuration in which it will be used.
- Refer to further guidance below regarding the EMC environment in which the device should be used.

Wireless communications equipment such as wireless home network devices, mobile phones, cordless telephones and their base stations, walkie-talkies can affect this equipment and should be kept at least a distance **d** away from the equipment. The distance **d** is calculated by the MANUFACTURER from the 800 MHz to 2,5 GHz column of Table 6 of IEC 60601-1-2,

**Recommended separation distance between portable and mobile RF communications equipment and the thermometer**

The thermometer are intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customers or the users of these thermometer can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the thermometer as recommended below, according to the maximum output power of the communications equipment.

Output Power of Transmitter in Watt	Separation distance according to frequency of transmitter in meter		
	150 kHz to 80 MHz $d = 1,2 \sqrt{P}$	80 MHz to 800 MHz $d = 1,2 \sqrt{P}$	800 MHz to 2,5 GHz $d = 2,3 \sqrt{P}$
0.01	0.12	0.12	0.23
0.1	0.38	0.38	0.78
1	1.2	1.2	2.3
10	3.8	3.8	7.3
100	12	12	23

For transmitters rated at a maximum output power not listed above, the recommended separation distance **d** in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where **P** is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

Note: At 80MHz and 800MHz, the separation distance for the higher frequency range applies

Note: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects, and people.

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© Miniland, S.A. 2011
Miniland S.A. P.Ind. La Marjal I C/ La Patronal 8-10
Shenzhen Dongdixin Technology Co., Ltd
Shanghai International Holding Corp. GmbH (Europe)